




DeskRay PT

- ✓ Fast & Easy to use
- ✓ Intuitive GUI design
- ✓ Reliable Application
- ✓ All Features Standard
- ✓ Budget Friendly
- ✓ Fully Integrated for TIMS



The **DeskRay PT** DUET offers the most essential features required by busy sites to produce universal exchange medical CDs & DVDs for referring physicians and patients alike. The **DeskRay PT** was designed with the end-user in mind, thereby creating a friendly and intuitive graphical user interface.



Each **DeskRay PT** DUET system provides true value and functionality. Standard features include the ability to merge patient studies, PACS Query & Retrieve, Burn CD/DVDs with spanning, DICOM Viewer of your choice, DICOM Image Importer.

DUET

The DUET 410 Recorder Printer uses a Plextor 760 CD-DVD combination drive with a recording speed up to 48x CD-R; 16x DVD±R. Utilizing the DUET ink cartridge, it prints up to 1200 DPI on Inkjet-printable CD-R/DVD±R media of both 8cm & 12cm sizes, with speeds ranging from 20 to 180 seconds.

only \$4,995



www.NautilusMedical.com
866-520-6477